

Wallaby Pediatric Wheelchair

HCPCS: E1236



FIGURE A

- Standard Transit Tie Down Brackets
- WC19 Compliant when used with WB8021 for users between 50 and 80 lbs
- Height adjustable push handles
- Flip back padded desk arms
- Pelvic belt standard
- Push to lock wheel locks
- 22" Rear wheels with 6" front casters
- Swingaway footrests with tool free height adjustment and calf strap

PRODUCTS	Description
WB1200-2GFR	Wallaby, 12", Fire Truck Red
WB1400-2GJB	Wallaby, 14", Jet Fighter Blue

SPECIFICATIONS	Description
Armrest Height	7"
Backrest Height (FC 2000N FC 4000N)	14"
Product Weight with footrests	37 lbs.
Seat To Floor Height	19"
Weight Capacity	150 lbs
Weight Capacity (13046 13046-HD)	150 lbs.
Weight Capacity when used in transport	50-80 lbs.

ACCESSORIES	Description
WB 8021	Headrest Extension, Wallaby, 1ea, EA
WB ELR	Elevating Leg Rests, Wallaby, 1pr, PR
STDS819	Anti Tippers, Cruiser, 1pr/bx, BX